


**New
Product...**

CovarisTM
the sample prep advantage

Covaris 96 microTUBETM Plate

Preassembled



Designed for Automation

| KEY FEATURES | BENEFITS |
|---|--|
| <ul style="list-style-type: none"> • Pre-assembled plate | <ul style="list-style-type: none"> ✓ Ready to use - no manual tube loading required |
| <ul style="list-style-type: none"> • Integrated, unitized plate assembly | <ul style="list-style-type: none"> ✓ Tubes permanently installed in plate ✓ Maintains uniform sample history and sample chain of custody |
| <ul style="list-style-type: none"> • No crimps, caps, or septa⁽¹⁾ | <ul style="list-style-type: none"> ✓ Low insertion/extraction force seal designed for automated pipetting |
| <ul style="list-style-type: none"> • Stackable, nested plates | <ul style="list-style-type: none"> ✓ Automation-ready |
| <ul style="list-style-type: none"> • Covaris microTUBE with AFA fiber | <ul style="list-style-type: none"> ✓ Uses existing E-series protocols |

(1) Replaces microTUBES with crimp caps (PN #520052) and 96 microTUBE rack (PN #500143)

Ordering Information

| Part Number | Product Name | Description |
|-------------|--------------------|--|
| 520051 | 96 microTUBE Plate | 6X16mm round bottom glass microTUBE in a 96 position microplate format |

Contact Information
Tel: +1 781-932-3959
Fax: +1 781-932-8705
Email: info@covarisinc.com
Web: www.covarisinc.com

USA Headquarters
Covaris, Inc.
14 Gill Street, Unit H
Woburn, Massachusetts 01801-1721
USA